



Glycoprotein

OVERVIEW

PRODUCT REFERENCE	PX-P4466
ORIGIN SPECIES	
FRAGMENT	Met1-Leu524
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	59.92kDa
FORM	
BUFFER	50 mM Tris-HCl pH 7.4, 150 mM NaCl, 10% glycerol, 4 mM EDTA, 1% CHAPS
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P08667
ALIAS/SYNONYMS	/

