



Undecaprenyl-diphosphatase(uppP)

OVERVIEW

PRODUCT REFERENCE	PX-P4520
ORIGIN SPECIES	
FRAGMENT	Gly70-Phe273
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	23.49kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P67389
ALIAS/SYNONYMS	Bacitracin resistance protein, Undecaprenyl pyrophosphate

