



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Up-regulated in infective sporozoites(PBANKA_0501200)

OVERVIEW

PRODUCT REFERENCE	PX-P4537
ORIGIN SPECIES	
FRAGMENT	Glu79-Ile220
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	18.09kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	A0A509AFU9
ALIAS/SYNONYMS	



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.