



Platelet-activating factor acetylhydrolase (PLA2G7)

OVERVIEW

PRODUCT REFERENCE	PX-P4546
ORIGIN SPECIES	
FRAGMENT	Met1~Asn441
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	50.94KDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q13093
ALIAS/SYNONYMS	PAFAH,1-alkyl-2-acetyl-glycerophosphocholine esterase,2-acetyl-1-alkyl-glycerophosphocholine esterase,gVIIA-PLA2,LDL-PLA(2),PAF 2-acylhydrolase

