



Alpha-enolase(Eno1)

OVERVIEW

PRODUCT REFERENCE	PX-P4549
ORIGIN SPECIES	
FRAGMENT	Met1~Lys434
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	48.8kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P06733
ALIAS/SYNONYMS	Eno-1,2-phospho-D-glycerate hydro-lyase,Enolase 1,NNE

