



Glycoprotein (G) gene from patient 07-83--GA-P

OVERVIEW

PRODUCT REFERENCE	PX-P4552
ORIGIN SPECIES	
FRAGMENT	Met1~Arg517
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	59.11kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q89015
ALIAS/SYNONYMS	

