



Protein A, von Willebrand factor binding protein Spa(spa)

OVERVIEW

PRODUCT REFERENCE	PX-P4556
ORIGIN SPECIES	
FRAGMENT	Met1~Lys327
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	38.33kDa
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	A0A098GZT2 or P38507
ALIAS/SYNONYMS	Immunoglobulin G-binding protein A, IgG-binding protein A, Staphylococcal protein A

