



Galectin-9(LGALS9)

OVERVIEW

PRODUCT REFERENCE	PX-P4585
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ala2-Thr323
TAG	N-terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	61.76kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O00182-2
ALIAS/SYNONYMS	LGALS9, Medium, Gal-9delta5, D5,Gal-9, Ecalectin, Tumor antigen HOM-HD-21

