



Serine-threonine-protein phosphatase 2A 65 kDa regulatory subunit A alpha isoform(Ppp2r1a)

OVERVIEW

PRODUCT REFERENCE	PX-P4598
ORIGIN SPECIES	
FRAGMENT	Met1-Ala589
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	66.4kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q76MZ3
ALIAS/SYNONYMS	PP2A subunit A isoform PR65-alpha, PP2A subunit A isoform R1-alpha

