



Transcription factor Sp5(SP5)

OVERVIEW

PRODUCT REFERENCE	PX-P4603
ORIGIN SPECIES	
FRAGMENT	Ile68-Ile274
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	24.5kDa
FORM	
BUFFER	PBS,pH7.5, 4Murea
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	T2MC68
ALIAS/SYNONYMS	Transcription factor Sp5

