



Epidermal Growth Factor proteins receptor(EGFR)

OVERVIEW

PRODUCT REFERENCE	PX-P4657
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gly696-Gly1022
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P00533
ALIAS/SYNONYMS	Proto-oncogene c-ErbB-1, Receptor tyrosine-protein kinase erbB-1, ERBB, ERBB1, HER1

