



Oncostatin-M-specific receptor subunit beta(OSMR)

OVERVIEW

PRODUCT REFERENCE	PX-P4681
ORIGIN SPECIES	
FRAGMENT	Met1-Met740
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	85.2kDa
FORM	Lyophilized
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q99650
ALIAS/SYNONYMS	Interleukin-31 receptor subunit beta,IL-31 receptor subunit beta,IL-31R subunit beta,IL-31R-beta,IL-31RB,OSMRB

