



Calbindin (CALB1)

OVERVIEW

PRODUCT REFERENCE	PX-P4697
ORIGIN SPECIES	
FRAGMENT	Met1-Asn261
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05937
ALIAS/SYNONYMS	Calbindin D28,D-28K,Vitamin D-dependent calcium-binding protein, avian-type,CAB27

