



Fibroblast Growth Factor proteins 21(FGF21)(29-209)

OVERVIEW

PRODUCT REFERENCE	PX-P4734
ORIGIN SPECIES	
FRAGMENT	His29-Ser209
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	20.61kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NSA1
ALIAS/SYNONYMS	FGF-21

