



Trans-acting T-cell-specific transcription factor GATA-3 (GATA3)

OVERVIEW

PRODUCT REFERENCE	PX-P4745
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Gly443
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	60.36kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P23771
ALIAS/SYNONYMS	GATA-binding factor 3,GATA3

