



Terminal nucleotidyltransferase 4B(Tent4b) Mouse

OVERVIEW

PRODUCT REFERENCE	PX-P4813
ORIGIN SPECIES	House Mouse
FRAGMENT	Phe2-Arg633
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	70.8kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q68ED3
ALIAS/SYNONYMS	Non-canonical poly(A) RNA polymerase PAPD5,Papd5,TRF4-2,TUTase 3,Topoisomerase-related function protein 4-2,Terminal uridylyltransferase 3,Terminal guanylyltransferase,PAP-associated domain-containing protein 5

