



Terminal nucleotidyltransferase 4A(TENT4A) Human

OVERVIEW

PRODUCT REFERENCE	PX-P4814
ORIGIN SPECIES	
FRAGMENT	Gln2-Arg772
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	83.4kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q5XG87
ALIAS/SYNONYMS	DNA polymerase sigma, LAK-1,

