



Fumarate hydratase, mitochondrial(FH)

OVERVIEW

PRODUCT REFERENCE	PX-P4826
ORIGIN SPECIES	
FRAGMENT	Ala45-Lys510
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	51.3kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P07954
ALIAS/SYNONYMS	Fumarase,HsFH

