



Fibronectin(FN1)AZN

OVERVIEW

PRODUCT REFERENCE	PX-P4840
ORIGIN SPECIES	
FRAGMENT	Glu1632-Ala1812
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	20.9kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P02751
ALIAS/SYNONYMS	Cold-insoluble globulin,CIG,FN

