



GAPDH protein - Glyceraldehyde-3-phosphate dehydrogenase

OVERVIEW

PRODUCT REFERENCE	PX-P4846
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Glu335
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P04406
ALIAS/SYNONYMS	GAPDH,Peptidyl-cysteine S-nitrosylase GAPDH,GAPD

