



CD40 ligand(CD40LG)

OVERVIEW

PRODUCT REFERENCE	PX-P4855
ORIGIN SPECIES	
FRAGMENT	Glu108-Leu261
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	18kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P29965
ALIAS/SYNONYMS	CD40-L,T-cell antigen Gp39,TNF-related activation protein,TRAP,Tumor necrosis factor ligand superfamily member 5,CD_antigen: CD154,CD40L, TNFSF5, TRAP

