



Interleukin-4 receptor subunit alpha(IL4R)

OVERVIEW

PRODUCT REFERENCE	PX-P4860
ORIGIN SPECIES	
FRAGMENT	Met1-His232
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	27.3kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P24394
ALIAS/SYNONYMS	IL4RA,IL-4 receptor subunit alpha,IL-4R subunit alpha,IL-4R-alpha,IL-4RA,CD_antigen:CD124

