



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Protein-arginine deiminase type-4(PADI4)

OVERVIEW

PRODUCT REFERENCE	PX-P4886
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1~Pro663
TAG	N-terminal Twin-Strep Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	77.29kDa
FORM	Liquid
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9UM07
ALIAS/SYNONYMS	PAD4, PADI5, PDI5, HL-60 PAD, Peptidylarginine deiminase IV, Protein-arginine deiminase type IV



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.