



## Recombinant Bovine Enterokinase(rb-EK)

### OVERVIEW

---

|                           |                    |
|---------------------------|--------------------|
| <b>PRODUCT REFERENCE</b>  | PX-P4889           |
| <b>ORIGIN SPECIES</b>     | Bovine             |
| <b>FRAGMENT</b>           |                    |
| <b>TAG</b>                | C-Terminal His Tag |
| <b>DELIVERY CONDITION</b> | Dry Ice            |

### PRODUCT INFORMATION

---

|                                |  |
|--------------------------------|--|
| <b>EXPRESSION SYSTEM</b>       | Prokaryotic expression                 |
| <b>HOST</b>                    | Escherichia coli (E.coli)              |
| <b>MOLECULAR WEIGHT</b>        | 26.2kDa                                |
| <b>FORM</b>                    | Frozen Liquid or Lyophilized           |
| <b>BUFFER</b>                  | 50 mM Tris-HCl, 50%Glycerol, pH 8.0    |
| <b>STABILITY &amp; STORAGE</b> | Store at -20 ° C, 12 months shelf life |
| <b>UNIPROT</b>                 |  |
| <b>ALIAS/SYNONYMS</b>          |  |

