



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Terminal Deoxynucleotidyl Transferase

OVERVIEW

PRODUCT REFERENCE PX-P4893
ORIGIN SPECIES
FRAGMENT
TAG
DELIVERY CONDITION Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM Prokaryotic expression
HOST Escherichia coli (E.coli)
MOLECULAR WEIGHT
FORM Frozen Liquid or Lyophilized
BUFFER 50mM K₂HPO₄, 100mM NaCl, 1mM DTT, 0.1% Tween20, 50% Glycerol, pH6,4
STABILITY & STORAGE Store at -20 ° C, 12 months shelf life
UNIPROT
ALIAS/SYNONYMS



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.