



Nuclease A (Benzonase)

OVERVIEW

PRODUCT REFERENCE	PX-P4894
ORIGIN SPECIES	<i>S. marcescens</i>
FRAGMENT	
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	<i>Escherichia coli</i> (E.coli)
MOLECULAR WEIGHT	29.9kDa
FORM	Frozen Liquid or Lyophilized
BUFFER	25 mM Tris-HCl pH 8.0, 20 mM NaCl, 2 mM MgCl ₂ , 50% glycerol
STABILITY & STORAGE	Store at -20 ° C, 12 months shelf life
UNIPROT	
ALIAS/SYNONYMS	NucA

