



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

IdeZ Protease

OVERVIEW

PRODUCT REFERENCE PX-P4895
ORIGIN SPECIES
FRAGMENT
TAG
DELIVERY CONDITION Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM Prokaryotic expression
HOST Escherichia coli (E.coli)
MOLECULAR WEIGHT 38 kDa
FORM Frozen Liquid or Lyophilized
BUFFER 50mM sodium phosphate, 150mM NaCl (Ph 6.6)
STABILITY & STORAGE Store at -20 ° C, 12 months shelf life
UNIPROT
ALIAS/SYNONYMS



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.