



Human CD42b-GP1BA recombinant protein, C-His Tag

OVERVIEW

PRODUCT REFERENCE	PX-P4904
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Leu531
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	58.52kDa
FORM	
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P07359
ALIAS/SYNONYMS	

