



Human CD66b recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5120
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gln35-Asp320
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	33.83 kDa
FORM	
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P31997
ALIAS/SYNONYMS	CD67 antigen, Carcinoembryonic antigen CGM6, Carcinoembryonic antigen-related cell adhesion molecule 8, CGM6, CD66b, Non-specific cross-reacting antigen NCA-95, CEACAM8

