



CD66a recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5131
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gln35-Gly428
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	43.34kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P13688
ALIAS/SYNONYMS	BGP, BGP1, CEACAM1

