



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

TGFB2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5188
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ala303-Ser414
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	12.32 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P61812
ALIAS/SYNONYMS	TGF-beta-2, Transforming Growth Factor proteins beta-2 proprotein



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.