



CD158j / KIR2DS2, N-GST, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5579
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	His22-Glu291
TAG	N-terminal GST Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	56.57 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P43631
ALIAS/SYNONYMS	

