



DSG3, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5686
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Confidential
TAG	No Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P32926
ALIAS/SYNONYMS	

