



IL27 Protein, C-His, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5787
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	(Phe 29-Pro 243) & (Met 1-Lys 229)
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	52 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8NEV9 & Q14213
ALIAS/SYNONYMS	

