



MYOM3 / Myomesin-3, N-His, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5852
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ser888-Glu1178
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	34.74 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q5VTT5
ALIAS/SYNONYMS	

