



OLFM4, N-GST & C-His, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5868
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Lys206-Gln510
TAG	N-terminal GST & C-His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	62.97 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q6UX06
ALIAS/SYNONYMS	

