



## Hydra vulgaris Transcription Factor Sp5-like, N-His, recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P5981
<b>ORIGIN SPECIES</b>	Hydra vulgaris
<b>FRAGMENT</b>	Ile45-Ile251(Leu56Ile, Glu70Asp, Ile100Val, Val221Ile)
<b>TAG</b>	N-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	24.55 kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	XP_004206818.1
<b>ALIAS/SYNONYMS</b>	

