



Monkeypox virus/MPXV A5L Recombinant Protein, C-His

OVERVIEW

PRODUCT REFERENCE	PX-P6078
ORIGIN SPECIES	Monkeypox virus (strain Zaire-96-I-16) (MPX)
FRAGMENT	Met1-Lys281
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	31.81 kDa
FORM	
BUFFER	Supplied as solution form in PBS pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	
ALIAS/SYNONYMS	A5L protein-like, A5L, Monkeypox virus (MPXV)

