



Human BMP2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6094
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gln283-Arg396
TAG	N-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	40.74 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P12643
ALIAS/SYNONYMS	BMP-2, BMP2A, Bone morphogenetic protein 2A, Bone morphogenetic protein 2, BMP2, BMP-2A

