



Nile tilapia NCCRP1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6115
ORIGIN SPECIES	Nile tilapia
FRAGMENT	Met1-Pro233
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	27.56 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	I3KFP2
ALIAS/SYNONYMS	NCCRP-1, NCCRP1, NCC receptor protein 1 homolog, F-box only protein 50, FBXO50, Non-specific cytotoxic cell receptor protein 1 homolog

