



HHV-3 -ORF68 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6118
ORIGIN SPECIES	HHV-3
FRAGMENT	Met1-Tyr538
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	61.64 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9J3M8
ALIAS/SYNONYMS	Envelope glycoprotein E, gE, RF68, Varicella zoster virus (VZV), HHV-3

