



Human CD100-SEMA4D recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6132
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Arg734
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	82.46 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	AAH54500.1
ALIAS/SYNONYMS	A8, Semaphorin-4D, GR3, BB18, SEMAJ, C9orf164, CD100, SEMA4D, Drug Target

