



Human IgG1 Fc region recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6139
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Asp104-Lys330
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	26 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01857
ALIAS/SYNONYMS	Ig gamma-1 chain C region EU, Immunoglobulin heavy constant gamma 1, IGHG1, Ig gamma-1 chain C region NIE, Ig gamma-1 chain C region, Ig gamma-1 chain C region KOL

