



Human CD85f-LILRA5 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6148
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gly42-Arg268
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	28.15 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	A6NI73
ALIAS/SYNONYMS	Leukocyte immunoglobulin-like receptor subfamily A member 5, Immunoglobulin-like transcript 11, ILT-11, Leukocyte immunoglobulin-like receptor 9, LIR9, CD85f, LIR-9, LILRB7, CD85 antigen-like family member F, LILRA5, ILT11

