



Human CD85a-LILRB3 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6158
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Glu443
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	49.49 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O75022
ALIAS/SYNONYMS	CD85a, ILT-5, Monocyte inhibitory receptor HL9, LILRB3, Immunoglobulin-like transcript 5, CD85 antigen-like family member A, ILT5, Leukocyte immunoglobulin-like receptor 3, LIR-3, Leukocyte immunoglobulin-like receptor subfamily B member 3, LIR3

