



Human CD85b-LILRA6 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6212
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Asn447
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	49.66 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q6PI73
ALIAS/SYNONYMS	Immunoglobulin-like transcript 8, ILT8, Leukocyte immunoglobulin-like receptor subfamily A member 6, LILRA6, Leukocyte Ig-like receptor, ILT-8

