



Romiplostim recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6248
ORIGIN SPECIES	
FRAGMENT	IEGPTLRQWLAARA
TAG	N-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	29.42 kDa
FORM	
BUFFER	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	
ALIAS/SYNONYMS	

