



Human ERBB4 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6273
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Cys29-Asn278
TAG	N-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	30.59 kDa
FORM	
BUFFER	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q15303
ALIAS/SYNONYMS	HER4, Tyrosine kinase-type cell surface receptor HER4, p180erbB4, s80HER4, 4ICD, Receptor tyrosine-protein kinase erbB-4, Proto-oncogene-like protein c-ErbB-4, E4ICD

