



4F2Mab Biosimilar - Anti-ganglioside GD2, SFRP1, SLC3A2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1073
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a-nd; IgG1-nd
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	4F2Mab, 14G2a, 4F2, ganglioside GD2, SFRP1, SLC3A2, anti-ganglioside GD2, SFRP1, SLC3A2

