



Afelimomab Biosimilar - Anti-TNFSF2, TNF-alpha, TNFA mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1088
CLONALITY	Monoclonal Antibody
ISOTYPE	F(ab') ₂ -G3-kappa
SOURCE	CAS 156227-98-4

PRODUCT INFORMATION

ORIGIN SPECIES	Mus musculus
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Afelimomab, MAK-195F, TNFSF2, TNF-alpha, TNFA, anti-TNFSF2, TNF-alpha, TNFA

